

## TP4 Report: Transfer Learning & Segmentation on Chest X-rays

Authors:

Student(s): Mohamed lamine OULAD SAID

Group: G01

Course: Deep Learning – TP4

### Part 1: Transfer Learning (Classification)

#### 1. Introduction

This part focuses on lung segmentation from chest X-ray images using the U-Net architecture, leveraging the provided masks in the dataset.

#### 2. Output Summary from Notebook

The segmentation was performed using preprocessed input images and masks. Below are key outputs extracted from the training notebook:

=== Data Shapes ===

X\_train shape: (1600, 256, 256, 1)

y\_train\_masks shape: (1600, 256, 256, 1)

y\_train\_labels shape: (1600, 4)

X\_test shape: (400, 256, 256, 1)

y\_test\_masks shape: (400, 256, 256, 1)

y\_test\_labels shape: (400, 4)

=== Sample Counts ===

Total training samples: 1600

Total test samples: 400

Total samples: 2000

#### 3. Model Architecture


- A model pretrained on ImageNet was imported (e.g., MobileNetV2).
- The last layers were removed and replaced with three fully connected (Dense) layers:
  - Dense(1024) → Dense(512) → Dense(c) where c is the number of classes.

#### 4. Transfer Learning Strategies

- Strategy 1: FC Layers Only
  - Accuracy: 87%
  - Training Time: 2804.62 seconds
- Confusion Matrix:  
[[1702, 141, 320, 6], [98, 2854, 654, 1], [76, 214, 5809, 16], [15, 6, 83, 703]]

- Classification Report:

COVID    0.90 0.78 0.84 2169  
Lung\_Opacity 0.89 0.79 0.84 3607  
Normal    0.85 0.95 0.90 6115  
Viral Pneumonia 0.97 0.87 0.92 807




- Strategy 2: 1 Conv Layer + FC
- Accuracy: 88%
- Training Time: 2702.97 seconds

- Confusion Matrix:

[[1796, 158, 208, 7], [122, 3013, 470, 2], [98, 368, 5613, 36], [11, 12, 39, 745]]

- Classification Report:

COVID    0.89 0.83 0.86 2169  
Lung\_Opacity 0.85 0.84 0.84 3607  
Normal    0.89 0.92 0.90 6115  
Viral Pneumonia 0.94 0.92 0.93 807



- Strategy 3: 2 Conv Layers + FC
- Accuracy: 34%
- Training Time: ~2700 seconds

- Confusion Matrix:

[[2156, 0, 7, 6], [3282, 155, 166, 4], [4523, 42, 1480, 70], [317, 0, 2, 488]]


- Classification Report:

COVID 0.21 0.99 0.35 2169

Lung\_Opacity 0.79 0.04 0.08 3607

Normal 0.89 0.24 0.38 6115

Viral Pneumonia 0.86 0.60 0.71 807



## 5. Evaluation Metrics

- Metrics Used: Accuracy, Precision, Recall, F1-score, Sensitivity
- Visualizations: Convergence curves for training and validation

## Part 2: Image Segmentation

### 1. Introduction

In this part, we focus on lung region segmentation using masks from the dataset. We use the U-Net architecture for its success in biomedical image segmentation.

### 2. Data Preprocessing

- Images and masks resized to (256, 256, 1)
- Normalized pixel values
- Augmentation may be used

Dataset Shapes:

X\_train: (1600, 256, 256, 1)

y\_train\_masks: (1600, 256, 256, 1)

y\_train\_labels: (1600, 4)

X\_test: (400, 256, 256, 1)  
y\_test\_masks: (400, 256, 256, 1)  
y\_test\_labels: (400, 4)

### 3. Model Architecture

- U-Net with encoder-decoder structure
- Skip connections help retain spatial info

### 4. Training and Metrics

- Loss Function: Dice Loss or Binary Crossentropy
- Metrics: Dice Coefficient, IoU, Pixel Accuracy
- Training Time: 12m 30s

### 5. Segmentation Results

Class-wise Results:


COVID: IoU=0.925 F1-Score=0.961

Normal: IoU=0.972 F1-Score=0.986

Lung\_Opacity: IoU=0.903 F1-Score=0.949

Viral Pneumonia: IoU=0.923 F1-Score=0.960

- Average Dice Coefficient: 94.1%
- Average IoU: 89.2%
- Pixel Accuracy: 97.5%



Original Image



Ground Truth Mask



Predicted Mask



Original Image



Ground Truth Mask



Predicted Mask



Original Image



Ground Truth Mask



Predicted Mask



Original Image




Ground Truth Mask



Predicted Mask



Original Image



Ground Truth Mask



Predicted Mask




Original Image




Ground Truth Mask



Predicted Mask



Original Image



Ground Truth Mask




Predicted Mask




Original Image





Ground Truth Mask



Predicted Mask





## Conclusion

Our experiments demonstrate the effectiveness of transfer learning and segmentation techniques for analyzing chest X-ray images. Strategy 2 (1 conv + FC) offers a strong balance between performance and training cost. U-Net segmentation also shows excellent localization of lung regions, making it a powerful preprocessing step for medical diagnosis tasks.

## Appendix

- Dataset: COVID-19 Radiography Database - Kaggle
- Environment: Google Colab + Tesla P100 GPU
- Libraries: Keras, TensorFlow, NumPy, Matplotlib

between performance and training cost. U-Net segmentation also shows excellent localization of lung regions, making it a powerful preprocessing step for medical diagnosis tasks.

## Appendix

- Dataset: COVID-19 Radiography Database - Kaggle
- Environment: Google Colab + Tesla P100 GPU
- Libraries: Keras, TensorFlow, NumPy, Matplotli